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1. Mehio R, Ruehle M, Catreux S, et al. DRAGEN Wins at Precision-FDA Truth Challenge V2 Showcase Accuracy Gains from Alt-aware Mapping and Graph Reference Genomes. [illumina.com/science/genomics-research/dragen-wins-precisionfda-challenge-showcase-accuracy-gains.html](https://illumina.com/science/genomics-research/dragen-wins-precisionfda-challenge-showcase-accuracy-gains.html). Published 2020. Accessed January 12, 2023.
2. Illumina. Accuracy improvements in germline small variant calling with the DRAGEN Bio-IT Platform. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/dragen-v4-accuracy-app-note-m-gl-01016/dragen-v4-accuracy-app-note-m-gl-01016.pdf](https://illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/dragen-v4-accuracy-app-note-m-gl-01016/dragen-v4-accuracy-app-note-m-gl-01016.pdf). Published 2022. Accessed February 1, 2023.
3. Pacific Biosciences. Preparing DNA for PacBio HiFi sequencing—Extraction and quality control. [pacb.com/literature/technical-note-preparing-dna-for-pacbio-hifi-sequencing-extraction-and-quality-control/](https://pacb.com/literature/technical-note-preparing-dna-for-pacbio-hifi-sequencing-extraction-and-quality-control/). Published 2022. Accessed January 12, 2023.
4. Pacific Biosciences. Preparing whole genome and metagenome libraries using SMRTbell prep kit 3.0. [pacb.com/wp-content/uploads/Procedure-checklist-Preparing-whole-genome-and-metagenome-libraries-using-SMRTbell-prep-kit-3.0.pdf](https://pacb.com/wp-content/uploads/Procedure-checklist-Preparing-whole-genome-and-metagenome-libraries-using-SMRTbell-prep-kit-3.0.pdf). Published 2022. Accessed October 5, 2023.
5. Oxford Nanopore Technologies. Ligation Sequencing Kit. <https://store.nanoporetech.com/us/productDetail/?id=ligation-sequencing-kit-v14>. Accessed October 5, 2023.
6. Pacific Biosciences. Low Yield Troubleshooting Guide. [pacb.com/wp-content/uploads/Guide-Low-Yield-Troubleshooting.pdf](https://pacb.com/wp-content/uploads/Guide-Low-Yield-Troubleshooting.pdf). Published 2018. Accessed January 12, 2023.
7. Illumina. Illumina Complete Long Read Prep, Human data sheet. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-data-sheet-m-gl-01420/illumina-long-read-prep-data-sheet-m-gl-01420.pdf](https://illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-data-sheet-m-gl-01420/illumina-long-read-prep-data-sheet-m-gl-01420.pdf). Published 2022. Accessed September 22, 2023.
8. Illumina. Comprehensive whole-genome sequencing with Illumina Complete Long Read Prep, Human technical note. [illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-tech-note-m-gl-01421/ilmn-long-read-hu-tech-note-m-gl-01421.pdf](https://illumina.com/content/dam/illumina/gcs/assembled-assets/marketing-literature/illumina-long-read-prep-human-tech-note-m-gl-01421/ilmn-long-read-hu-tech-note-m-gl-01421.pdf). Published 2022. Accessed September 22, 2023.
9. Bekritsky MA, Colombo C, Eberle MA. Identifying genomic regions with high quality single nucleotide variant calling. [illumina.com/science/genomics-research/articles/identifying-genomic-regions-with-high-quality-single-nucleotide-.html](https://illumina.com/science/genomics-research/articles/identifying-genomic-regions-with-high-quality-single-nucleotide-.html). Published 2021. Accessed August 30, 2023.
10. Wagner J, Olson ND, Harris L, et al. Curated variation benchmarks for challenging medically relevant autosomal genes. *Nat Biotechnol.* 2022;40(5):672–680. doi:10.1038/s41587-021-01158-1
11. PharmGKB. VIPs: Very Important Pharmacogenes. [pharmgkb.org/vips](https://pharmgkb.org/vips). Accessed September 22, 2023.
12. National Library of Medicine. GTR: Genetic Testing Registry. Precision HealthPGx Panel (25 Genes). [ncbi.nlm.nih.gov/gtr/tests/593428/](https://ncbi.nlm.nih.gov/gtr/tests/593428/). Updated November 29, 2022. Accessed September 22, 2023.
13. Pratt VM, Everts RE, Aggarwal P, et al. Characterization of 137 Genomic DNA Reference Materials for 28 Pharmacogenetic Genes: A GeT-RM Collaborative Project. *J Mol Diagn.* 2016;18(1):109–123. doi:10.1016/j.jmoldx.2015.08.005
14. Miller DT, Lee K, Abul-Husn NS, et al. ACMG SF v3.1 list for reporting of secondary findings in clinical exome and genome sequencing: A policy statement of the American College of Medical Genetics and Genomics (ACMG). *Genet Med.* 2022;24(7):1407–1414. doi:10.1016/j.gim.2022.04.006
15. Kulski JK, Suzuki S, Shiina T. Human leukocyte antigen super-locus: nexus of genomic supergenes, SNPs, indels, transcripts, and haplotypes. *Hum Genome Var.* 2022;9(1):49. doi:10.1038/s41439-022-00226-5
16. Roessler K. Illumina Complete Long Reads software analysis workflow for human WGS. [illumina.com/science/genomics-research/articles/complete-long-read-software-analysis.html](https://illumina.com/science/genomics-research/articles/complete-long-read-software-analysis.html). Published 2023. Accessed September 22, 2023.

| 製品情報                                                                                     |          | 製品情報                                                                     |          |
|------------------------------------------------------------------------------------------|----------|--------------------------------------------------------------------------|----------|
| 製品                                                                                       | カタログ番号   | 製品                                                                       | カタログ番号   |
| Illumina Complete Long Read Prep with Enrichment, Human (24 samples)                     | 20113832 | Illumina Unique Dual Indexes, LT (48 indexes, 48 samples)                | 20098166 |
| Illumina Complete Long Read Prep with Enrichment, Human (96 samples)                     | 20113833 | Illumina DNA/RNA UD Indexes Set A, Tagmentation (96 indexes, 96 samples) | 20091654 |
| Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel (24 samples) | 20113834 | Illumina DNA/RNA UD Indexes Set B, Tagmentation (96 indexes, 96 samples) | 20091656 |
| Illumina Complete Long Read Prep with Enrichment, Human Comprehensive Panel (96 samples) | 20113835 | Illumina DNA/RNA UD Indexes Set C, Tagmentation (96 indexes, 96 samples) | 20091658 |
| Illumina Human Comprehensive Panel (24 samples)                                          | 20113836 | Illumina DNA/RNA UD Indexes Set D, Tagmentation (96 indexes, 96 samples) | 20091660 |
| Illumina Human Comprehensive Panel (96 samples)                                          | 20113837 | Illumina Analytics - 1 iCredit                                           | 20042038 |
| Illumina Custom Enrichment Panel v2 (32 µl, 120 bp)                                      | 20073953 |                                                                          |          |
| Illumina Custom Enrichment Panel v2 (384 µl, 120 bp)                                     | 20073952 |                                                                          |          |
| Illumina Custom Enrichment Panel v2 (1536 µl, 120 bp)                                    | 20111339 |                                                                          |          |

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© 2024 Illumina, Inc. All rights reserved.